



## QV Cream - 500g (tub)

**\$20.15**

In stock  
Only **500** left  
**SKU** 11062



 Add to Wish List

 Add to Compare

Highly concentrated moisturising cream for dry skin.



### Details

This rich, protective moisturising cream replenishes the skin and helps to hydrate dry areas. It contains squalene, a natural component found in the skin's oil. Areas prone to dryness like the elbows, knees and hands will just love what QV Cream can do to help.

### Benefits:

- Recommended for moisturising dry skin.
- Helps to seal in natural moisture by acting as a protective barrier.
- Helps protect against dryness.
- pH balanced, low-irritant formulation, free from fragrance, colour, lanolin and propylene glycol.
- Suitable for skin conditions such as dry, scaly skin, eczema, psoriasis and dermatitis.
- Suitable for sensitive skin.
- Suitable for infants and the elderly.
- Vegan Friendly - contains no animal derived ingredients.



## Directions:

Apply to affected areas as often as required, particularly after bathing, showering, shaving, exposure to harsh climatic conditions and at night.

## Cautions:

void eyes. For external use only. Store below 30°C.

## Ingredients:

Aqua (Water), Paraffinum Liquidum, Glycerin, Petrolatum, Cetearyl Alcohol, Dimethicone, Squalane, Glyceryl Stearate SE, Cetareth-20, Stearic Acid, Laureth-3, Glyceryl Stearate, Methylparaben, Dichlorobenzyl Alcohol.