



QV Cream - 500g (pump)

\$21.60

In stock
Only **500** left
SKU 10123

 Add to Wish List

 Add to Compare

Highly concentrated moisturising cream for dry skin.



Details

This rich, protective moisturising cream replenishes the skin and helps to hydrate dry areas. It contains squalene, a natural component found in the skin's oil. Areas prone to dryness like the elbows, knees and hands will just love what QV Cream can do to help.

Benefits:

- Recommended for moisturising dry skin.
- Helps to seal in natural moisture by acting as a protective barrier.
- Helps protect against dryness.
- pH balanced, low-irritant formulation, free from fragrance, colour, lanolin and propylene glycol.
- Suitable for skin conditions such as dry, scaly skin, eczema, psoriasis and dermatitis.
- Suitable for sensitive skin.
- Suitable for infants and the elderly.
- Vegan Friendly - contains no animal derived ingredients.



Directions:

Apply to affected areas as often as required, particularly after bathing, showering, shaving, exposure to harsh climatic conditions and at night.

Cautions:

void eyes. For external use only. Store below 30°C.

Ingredients:

Aqua (Water), Paraffinum Liquidum, Glycerin, Petrolatum, Cetearyl Alcohol, Dimethicone, Squalane, Glyceryl Stearate SE, Cetareth-20, Stearic Acid, Laureth-3, Glyceryl Stearate, Methylparaben, Dichlorobenzyl Alcohol.