



AEROWOUND Combine Dressing 10 x 20cm Bag/50

\$72.46

In stock
SKU ACD1020S

 Add to Wish List

 Add to Compare

Combine Dressing bag/50

Details

Size: 10 x 20cm **Pricing Unit:** bag/50 **Carton QTY:** 8 The AeroWound Combine Dressing is recommended where a highly absorbent dressing is required for large and heavily bleeding wounds.

Low Adherent

- The AeroWound Combine Dressing has a special low adherent coating which minimizes trauma associated with changing wound dressings

Hypo-Allergenic

- The AeroWound Combine Dressing is hypo-allergenic with a low risk of wound irritation

Highly Absorbent

- The AeroWound Combine Dressing is capable of absorbing fluids from large and heavily bleeding wounds to stop bleeding and reduce infection risk.

Individually packed and sterile

- Dressing packs are available individually on demand in an easy peel sterile pouch to ensure wound cleanliness.



Range of Sizes Available

- Ranges of dressing sizes are available to suit both small and large sized wounds effectively.





More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes