



## AEROWOUND Combine Dressing 20 x 40cm

**\$8.66**

In stock  
SKU ACD20405

 Add to Wish List

 Add to Compare

Combine Dressing dressing/1



### Details

**Size:** 20 x 40cm **Pricing Unit:** unit/1 **Carton QTY:** 100 The AeroWound Combine Dressing is recommended where a highly absorbent dressing is required for large and heavily bleeding wounds.

### Low Adherent

- The AeroWound Combine Dressing has a special low adherent coating which minimizes trauma associated with changing wound dressings

### Hypo-Allergenic

- The AeroWound Combine Dressing is hypo-allergenic with a low risk of wound irritation

### Highly Absorbent

- The AeroWound Combine Dressing is capable of absorbing fluids from large and heavily bleeding wounds to stop bleeding and reduce infection risk.

### Individually packed and sterile

- Dressing packs are available individually on demand in an easy peel sterile pouch to ensure wound cleanliness.



## Range of Sizes Available

- Ranges of dressing sizes are available to suit both small and large sized wounds effectively.





### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes