



AEROCREPE Light Cotton Crepe Bandage 7.5cm x 4M Wrap/12

\$46.24

In stock
SKU ACL75

 Add to Wish List

 Add to Compare

Light Crepe Bandage 7.5cm


Details

Size: 7.5cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 30 The AEROCREPE Bandage is a high quality cotton crepe bandage, recommended for low support requirements and can be used for the retention of dressings.

Light Crepe

AEROCREPE light cotton crepe bandages are of a standard conforming nature which is durable and suitable for dressing retention and light support requirements.

Easy to Open Wrapper

AEROCREPE light cotton crepe bandages have a easy open tab built into the wrapper of the bandage. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes