



AEROCREPE Medium Cotton Crepe Bandage 15cm x 4M Wrap/12

\$87.08

In stock
SKU ACM150

 Add to Wish List

 Add to Compare

Medium Crepe Bandage 15cm



Details

Size: 15cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 20 The AEROCREPE Bandage is a high quality cotton crepe bandage, recommended for low support requirements and can be used for the retention of dressings.

Medium Crepe

AEROCREPE medium cotton crepe bandages are of a standard conforming nature which is durable and suitable for dressing retention and medium support requirements.

Easy to Open Wrapper

AEROCREPE medium cotton crepe bandages have a easy open tab built into the wrapper of the bandage. 



More Information

| | |
|------------------------|----------|
| MANUFACTURER_MS | Germany |
| MODEL NUMBER_MS | K61SD521 |
| ECOLABEL_MS | yes |