



## AEROCREPE Medium Cotton Crepe Bandage 5cm x 4M Wrap/12

**\$36.92**

In stock  
**SKU** ACM50

 Add to Wish List

 Add to Compare

Medium Crepe Bandage 5cm

### Details

**Size:** 5cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 40 The AEROCREPE Bandage is a high quality cotton crepe bandage, recommended for low support requirements and can be used for the retention of dressings.

### Medium Crepe

AEROCREPE medium cotton crepe bandages are of a standard conforming nature which is durable and suitable for dressing retention and medium support requirements.

### Easy to Open Wrapper

AEROCREPE medium cotton crepe bandages have a easy open tab built into the wrapper of the bandage. 



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes