




Small Stainless Steel Trolley with Drawer 50 x 50 x 97cm

\$3,265.56

In stock
SKU ACR12

 Add to Wish List

 Add to Compare
Trolley with Drawer

Details

Size: 50 x 50 x 97cm **Pricing Unit:** unit/1 **Carton QTY:** 1

Medical/Instrument Trolleys

These units provide a sturdy surface and storage for medical supplies of instruments. Stainless steel construction makes for a durable, aesthetic and easy to disinfect finish. Supplied flat packed with minor assembly required.

Product Features:

- 1 Drawer: 125mm deep
- High quality 304 grade stainless steel
- Safety fold to all edges
- Rubber wheels, front 2 locking
- Full extension telescopic ball bearing drawer runners



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes