



# AERODIAGNOSTIC Pistol Grip Infrared Forehead Thermometer

**\$144.25**

In stock  
**SKU** ADT20

 Add to Wish List

 Add to Compare

AERODIAGNOSTIC Pistol Grip Infrared Forehead Thermometer

## Details

**Colour:** Blue **Pricing Unit:** unit/1 **Carton QTY:** 40

ARTG No.334356

The AERODIAGNOSTIC range offers a selection of commonly used diagnostic supplies in the healthcare industry.

## Infrared Thermometers

AERODiagnostic Clinical Forehead Infrared Thermometers ensure accurate and stable body temperature readings, without the need for physical contact from a distance of around 5cm.

Simply point the thermometer towards the forehead to read the temperature shown on the LCD screen.

The AERODiagnostic Non-Contact Forehead IR Thermometer is specially designed to take the body temperature of a person regardless of room temperature.

## Specifications:

Product name: Non-contact Thermometer Battery: DC 3V (DC 1.5V AAA 2PCS) Size: 148mm (L)X45mm(W)X90mm(H)  
Weight: About 85g (with the battery) Accuracy: Body mode: 0.2(35.0~42.0) Object mode: 0.5(0~100) Measuring range:  
Body mode: 32.0~ 42.9 Object mode:0~100 Storage/transportation temperature: -20 to 55 Relative humidity<85%RH  
Measuring distant: 3-5cm



## More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes