



AERODIAGNOSTIC Pistol Grip Infrared Forehead Thermometer with Carry Case

\$280.58

In stock
SKU ADT25

 Add to Wish List

 Add to Compare

AERODIAGNOSTIC Pistol Grip Infrared Forehead Thermometer with Carry Case

Details

Colour: Purple **Pricing Unit:** unit/1 **Carton QTY:** 40 ARTG No.333667 The AERODIAGNOSTIC range offers a selection of commonly used diagnostic supplies in the healthcare industry.

Infrared Thermometers

AERODiagnostic Clinical Forehead Infrared Thermometers ensure accurate and stable body temperature readings, without the need for physical contact from a distance of around 5cm.

Simply point the thermometer towards the forehead to read the temperature shown on the LCD screen.

The AERODiagnostic Non-Contact Forehead IR Thermometer is specially designed to take the body temperature of a person regardless of room temperature.

Specifications:

Product name: Non-contact Thermometer Measuring Range: >35.0C to 42.0C / 95.0F to 107.6F Resolution: 0.1C /0.1F Basic Accuracy: 35.0C to 42.0C : 0.2C Measuring Distance: 5cm Response Time: 1 Second Size: 95mmL x 42mmW x 139mmH Weight: 85g (without batteries) Warranty: 12 Months Memorization: 32 SETS Battery: DC3V (AA*2) Relative Humidity: 85% Environment Temp: 16C to 35C



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes