



AEROFORM Conforming Bandage 5cm x 4M Wrap/12

\$15.52

In stock
SKU AF50

 Add to Wish List

 Add to Compare

Conforming Bandage 5cm



Details

Size: 5cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 60 The AeroForm Conforming bandage is a quality cost-effective bandage with excellent conforming attributes. Our Heavyweight and Premium range offers greater support and compression and are recommended for sprains, strains and snake bites.

Conforming Specifications

AeroForm Conforming Bandages are of a standard conforming nature which is durable and suitable for dressing retention.

Easy to Open Wrapper

AeroForm Conforming Bandages have a easy open tab built into the wrapper of the bandage. 



More Information

| | |
|------------------------|----------|
| MANUFACTURER_MS | Germany |
| MODEL NUMBER_MS | K61SD521 |
| ECOLABEL_MS | yes |