



## RESPONDER 4 Series Plastic Tacklebox First Aid Kit 42 x 21 x 22cm

**\$271.76**

In stock  
**SKU** AFAK4T



 Add to Wish List

 Add to Compare

RESPONDER 4 Series Plastic Tacklebox First Aid Kit 42 x 21 x 22cm



### Details


**Case:** Plastic

**Size:** 42 x 21 x 22cm

**Pricing Unit:** unit/1

The Responder Tackle Box Kit ensures workplace compliance for first aid, containing a wide scope of products to cater for any workplace emergencies. Internal compartments and shelves keep products organised and easy to find, assisting with prompt response.

### Available in 5 case options

Respond promptly to incidents with a variety of outer cases to suit all types of environments. This kit is available in 5 case options: a [Robust](#) rugged plastic case, Versatile [Softpack](#) reinforced fabric bag, a [Weather resistant](#) plastic case and a [Neat](#) plastic case - all our cases and bags are available with custom labelling. The softpack style is ideal for mobile workers, while the durable plastic cases are recommended for industrial areas, and all plastic and metal case options are able to be wall mounted. 

### Instant Restock

This First Aid Kit comes fully stocked and as part of our **Instant Restock™** system you can refill and service your client's kits quickly with refill pouches designed to fit into each case type. They provide quick access to all products through a



unique centre slit pocket without the need of unpacking each individual product.

### Contents meet Australian Standard AS2675-1983

#### KIT CONTENTS:

CODE	DESCRIPTION	QUANTITY
AP10125-I	AEROPLAST Plastic Plaster 72x19mm	25 pc
AFP50125-I	AEROPLAST Fabric Premium Strip 72mm x 19mm	25 pc
AB3-I	AEROBURN Burn Gel Sachet 3.5g	5 pc
AAS50-I	AEROAID Antiseptic Spray 50ml	1 pc
AFS001-I	AEROSHIELD CPR Face Shield	1 pc
AGNPF02-I	AEROGLOVE Nitrile Examination Gloves	5 Pairs
AGS753S-I	AEROSWAB Gauze Swabs 7.5cm x 7.5cm 3's	5 Pks
AW1000-I	AEROWASH Eye Wash Ampoule 15ml	8 pc
AW8001-I	AEROWIPE Cleansing Wipe	10 pc
ASP001-I	AEROPROBE Splinter Probes 3.7cm	10 pc
APD50S-I	AEROPAD Low Adherent Dressing 5cm x 5cm	6 pc
APD100S-I	AEROPAD Low Adherent Dressing 7.5cm x 10cm	3 pc
APD101S-I	AEROPAD Low Adherent Dressing 10cm x 10cm	2 pc
AF50-I	AEROFORM Conforming Bandage 5cm x 4M	3 Roll
AF75-I	AEROFORM Conforming Bandage 7.5cm x 4M	3 Roll
AFHP100-I	AEROFORM Heavyweight Conforming Bandage 10cm x 4m	1 Roll
AP325-I	AEROPORE Microporous Tape 2.5cm x 5M	1 Roll
AWD14S-I	AEROWOUND BPC Wound Dressing #14	1 Roll
AWD15S-I	AEROWOUND BPC Wound Dressing #15	1 Roll
ACD1020S-I	AEROWOUND Combine Dressing 10cm x 20cm	1 pc
ABC10-I	AEROBAND Triangular Bandage 110 x 110 x 155cm	2 pc
ATB130-I	AERORESCUE Emergency Rescue Blanket	1 pc
AEP1S-I	AEROPAD Eye Pads	4 pc
AII2500-I	AEROCOOL Instant Ice Pack 80g	1 pc
ASC13-I	AEROINSTRUMENTS Scissor 13cm	1 pc
AF13-I	AEROINSTRUMENTS Tweezer 11cm	1 pc
AAB100-I	AEROSUPPLIES Amputated Parts Bags Assorted Sizes	3 pc
ASN096P-I	AEROSUPPLIES Notebook & Pen	1 pc
ASP12-I	AEROPINS Safety Pin	12 pc
AGL001-I	AEROGUIDE First Aid Leaflet	1 pc



## More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes