



AEROINSTRUMENTS Stainless Steel Magnifying Glass Forceps 11cm

\$23.80

In stock
SKU AFG11

 [Add to Wish List](#)

 [Add to Compare](#)

Forceps S/S with Moving Head Glass

Details

Size: 11cm **Pricing Unit:** unit/1 **Carton QTY:** 750 The AEROINSTRUMENT range offers a wide range of surgical instruments manufactured from high quality stainless steel components to suit first aid requirements. Assorted sizes and specifications have been selected for use in common first aid procedures.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes