





AEROFILM PLUS Waterproof Island Dressing 10 x 12cm Box/3 (GST FREE)

\$7.99

In stock
SKU AFPW0103

 Add to Wish List
 Add to Compare
Film Dressing Pad

Details

Size: 10 x 12cm **Pricing Unit:** env/3 **Carton QTY:** 96 This post-operative dressing combines the advantages of AEROFILM with a SupaSorb™ absorbent and low adherent dressing pad to provide a superior dressing for post operative and delicate wound sites.

SupaBond™ Hypo-allergenic Adhesive

- Exclusively engineered for unrivaled adhesion, this hypoallergenic adhesive has a low irritation risk while ensuring the dressing stays in place over body contours.

Transparent AeroFilm

- The transparent AeroFilm allows the user to monitor the level of fluid and check for clinical signs of infection without removing the dressing. This allows for maximum usage of the dressing before removal.

Simple Application

- AeroFilm has an integrated frame applicator that provides a simple and effective application system. The frame enables the dressing to be applied crease-free and ensures adhesion. AeroFilm transparency enables correct positioning over the wound site.



Contraindications / Precautions

- AeroFilm should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes