



AEROFORM Size A Newborn Elastic Tubular Bandage 4cm x 10M

\$45.41

In stock
SKU AGA10

 [Add to Wish List](#)

 [Add to Compare](#)

Size A Tubular Bandage 10M




Details

Size: 4cm x 10M **Pricing Unit:** unit/1 **Carton QTY:** 1 AEROFORM Tubular Bandages are an elasticised support bandage suited to treatment of sports injuries, dislocations and sprains. The Size A tubular bandage is suitable for use on infant limbs.

Self Supporting

AEROFORM Tubular Bandages are held in place without tying or tapes, resulting in a comfortable long wearing bandage.

Provides Even Support

The AEROFORM Tubular Bandages provides even support across damaged joints and muscles due to its elastic self supporting properties resulting in a comfortable support bandage. Available in 1m lengths or 10m rolls. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes