





## AEROSWAB Non-Sterile White Gauze Swab 7.5 x 7.5cm Pack/100 (8PLY)

**\$4.26**

In stock  
**SKU** AGS75

-  Add to Wish List
  -  Add to Compare
- 8 ply Gauze Swab



### Details

**Size:** 7.5 x 7.5cm **Pricing Unit:** pack/100 **Carton QTY:** 80 The AEROSWAB™ Gauze Swab is highly absorbent cotton gauze swab that when soaked in fluid becomes an ideal item to clean and prepare wounds for dressing.

### Cotton Gauze

- AEROSWAB™ cotton gauze has a balanced combination of softness, abrasive texture and absorbency for effectively cleaning wounds.

### Multiple Sizes Available

- AEROSWAB™ gauze swabs are available in a number of different sizes to suit any application





## More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes