





## AEROSWAB Sterile White Gauze Swab 7.5 x 7.5cm Box/25 (Packs of 3)

**\$7.22**

In stock  
**SKU** AGS753S

-  Add to Wish List
  -  Add to Compare
- Gauze Swab - 3 Pcs

### Details

**Size:** 7.5 x 7.5cm **Pricing Unit:** box/25 **Carton QTY:** 24 The AEROSWAB™ Gauze Swab is highly absorbent cotton gauze swab that when soaked in fluid becomes an ideal item to clean and prepare wounds for dressing.

### Cotton Gauze

- AEROSWAB™ cotton gauze has a balanced combination of softness, abrasive texture and absorbency for effectively cleaning wounds.

### Multiple Sizes Available

- AEROSWAB™ gauze swabs are available in a number of different sizes to suit any application.

### Sterile

- AEROSWAB™ gauze swabs are also available in a sterile form for minimising infection risk when dressing clean wounds





### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes