



AEROCASE Small/Medium Metal Cabinet 29 x 42 x 17cm

\$96.41


In stock
SKU AKM002

 Add to Wish List

 Add to Compare

Metal Cabinet Small/Med Side Opening

Details

Size: 42 x 29 x 17cm **Pricing Unit:** cabinet **Carton QTY:** 1 The AeroCase range of wall mountable cabinets cater to a broad range of first aid kit requirements in the industry, and includes both metal and plastic constructions with multiple sizes and configurations. Cabinet size should reflect the size and capacity requirement for the specific workplace or industrial kit. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes