




## AEROCASE Medium Metal Cabinet 35 x 45 x 17cm

**\$112.31**


In stock  
**SKU** AKM003

 Add to Wish List

 Add to Compare

Metal Cabinet Medium Side Opening

### Details

**Size:** 35 x 45 x 17cm **Pricing Unit:** cabinet **Carton QTY:** 1 The AeroCase range of wall mountable cabinets cater to a broad range of first aid kit requirements in the industry, and includes both metal and plastic constructions with multiple sizes and configurations. Cabinet size should reflect the size and capacity requirement for the specific workplace or industrial kit. 



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes