



AEROCASE Medium Drop Front Metal Cabinet 43 x 31 x 15cm

\$92.41

In stock
SKU AKM006


 [Add to Wish List](#)

 [Add to Compare](#)

Metal Cabinet Medium Drop Front



Details

Size: 31 x 43 x 15cm **Pricing Unit:** cabinet **Carton QTY:** 1 The AeroCase range of wall mountable cabinets cater to a broad range of first aid kit requirements in the industry, and includes both metal and plastic constructions with multiple sizes and configurations. Cabinet size should reflect the size and capacity requirement for the specific workplace or industrial kit. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes