



AEROCASE Small White Plastic Cabinet with Key Latch 26 x 32 x 14cm

\$61.20

In stock
SKU AKP0031

 [Add to Wish List](#)

 [Add to Compare](#)
Plastic Cabinet White



Details

Size: 26 x 32 x 14cm **Pricing Unit:** cabinet **Carton QTY:** 1 The AeroCase range of wall mountable cabinets cater to a broad range of first aid kit requirements in the industry, and includes both metal and plastic constructions with multiple sizes and configurations. Cabinet size should reflect the size and capacity requirement for the specific workplace or industrial kit. This cabinet has a key latch.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes