



AEROTAPE White Microporous Paper Tape 1.25cm x 5M Box/24

\$13.73

In stock
SKU AP312



 Add to Wish List

 Add to Compare

Microporous Paper Tape



Details

Size: 1.25cm x 5M **Pricing Unit:** box/24 **Carton QTY:** 15 AEROPORE Microporous Paper Tape is a flexible and breathable paper tape suitable for securing dressings, especially to facial and fragile skin. AeroPore Non-woven Microporous Paper Tape is recommended for the securing of dressings and medical equipment on sensitive and fragile skin including facial areas.

Hypoallergenic

- AEROPORE Microporous Paper Tape has a very low irritant risk, and is suitable for use on fragile and delicate skin.

Individually Wrapped

- AEROPORE Microporous Paper Tape is available in boxes of individually wrapped rolls ensuring the tape is kept in optimum condition.





More Information

| | |
|------------------------|----------|
| MANUFACTURER_MS | Germany |
| MODEL NUMBER_MS | K61SD521 |
| ECOLABEL_MS | yes |