



AEROTAPE White Microporous Paper Tape with Dispenser 2.5cm x 9.1M Box/12

\$37.33

In stock
SKU AP325D

 Add to Wish List

 Add to Compare

Microporous Paper Tape with Dispenser



Details

Size: 2.5cm x 9M **Pricing Unit:** box/12 **Carton QTY:** 30 AEROPORE Microporous Paper Tape is a flexible and breathable paper tape suitable for securing dressings, especially to facial and fragile skin. AeroPore Non-woven Microporous Paper Tape is recommended for the securing of dressings and medical equipment on sensitive and fragile skin including facial areas.

Hypoallergenic

- AEROPORE Microporous Paper Tape has a very low irritant risk, and is suitable for use on fragile and delicate skin.

Ease of Use

- AEROPORE Microporous Paper Tape is available in a ready to use dispenser for easy dispensing and tear off

 Includes Tape Dispenser.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes