



## AEROPLAST Plastic X-Wide Strip 7.5 x 2.5cm Box/100

**\$7.07**

In stock  
**SKU AP801**

 Add to Wish List

 Add to Compare

AEROPLAST Extra Wide Strip



### Details

**Size:** 7.2 x 2.5cm **Pricing Unit:** box/100 **Carton QTY:** 36 The AEROPLAST Plastic dressing is a wash proof sheer bandage recommended for all general applications.

### SupaSorb™ Wound Pad

- An embossed wound pad provides superior wound management by draining blood or exudates into pockets in the wound pad. The Super absorbent wound pad minimises the number of dressing changes whilst protecting the wound from further injury. SupaSorb™ leaves the wound site clean and comfortable, while the non-adherent film reduces the pain and site trauma when changing the dressing.

### SupaBond™ Adhesive


- SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling.

### Conforming and Air Permeable

- The AeroPlast Plastic Dressing is made of unique materials that bend, stretch and conform to different shapes and forms, as well as being air permeable which allows the skin and wound to breathe as it heals.



### Contraindications / Precautions

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician. **This product is not a substitute for direct medical attention from a healthcare practitioner.** 



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes