



## AEROPLAST Rigid Sports Tape 3.8cm x 13.7M Box/12

**\$70.00**

In stock  
**SKU APRT38**



 Add to Wish List

 Add to Compare

Rigid Sports Tape Medium



### Details

**Size:** 3.8cm x 13.7M **Pricing Unit:** box/12 **Carton QTY:** 6

### Our New Premium Sports Tape Range

- Supabond™Hypoallergenic Adhesive
- High-Quality Performance Based Materials
- Competitive Price
- Latex Free
- Lightweight for Comfort and Performance
- Flexible Kinesiology and Rigid Options
- Trusted on and off the field

**AeroPlast Kinesiology Tapes** are a highly flexible elastic tape that adheres tightly to skin promoting increased blood flow and relieves pressure on injured muscles. Supabond™Hypoallergenic Adhesive is latex free with a low allergy risk that provides a secure adhesive with excellent conformability for joints in awkward areas.

**AeroPlast Rigid Sports Tapes** are a lightweight and latex free rigid tape, ideal for strapping and supporting strains and sprains. The tape is designed to be hand tearable.





### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes