



Green Terylene/Cotton Carry Sheet



\$79.61

In stock
SKU AR1017

 Add to Wish List

 Add to Compare

Green Cotton Carry Sheet Terylene



Details

Size: 190 x 68cm **Pricing Unit:** unit/1 **Carton QTY:** 1 The green cotton terylene carry sheet is a compact and convenient carry sheet suitable for a range of actions for fast patient recovery. This product is highly suitable for public emergency with its key attributes being flexibility and strength. The carry sheet folds to a small package; yet has 8 strong carry handles for patient handling.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes