



AEROSPLINT Junior Aluminium Foam Splint 45 x 11cm

\$10.52

In stock
SKU AR1045

 Add to Wish List

 Add to Compare

AEROSPLINT Junior Aluminium Foam Splint 45 x 11cm

Details

Size: 45 x 11cm (18" x 4.5") **Pricing Unit:** unit/1 **Carton QTY:** 100 The Junior aluminium foam splint is ideal for children because of the shorter length in comparison to the AR1018F/AR1018R adult foam splints.

Aluminium Foam Splints

An aluminium foam splint is a highly adaptive splint that can be rolled folded or formed into almost any application with the use of a quality support bandage.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes