




AEROWOUND Finger Dressing Env/3

\$4.21

In stock
SKU AWF03

 Add to Wish List

 Add to Compare
Finger Dressing

Details

Size: 4.5 x 4.5cm **Pricing Unit:** Env/3 **Carton QTY:** 96 The AEROWOUND Wound Dressings are a general purpose wound dressing for use in preventing infection and controlling bleeding. Available in sizes to suit wounds from finger injuries to large wounds.

Unique Finger Dressing

- The AEROWOUND Finger dressing is unique in the industry for its ability to treat finger wounds too large to cover with an adhesive dressing. This unique product combines a 4.5cm square low adherent and absorbent pad, with a wrapping leader and secure adhesive tab to protect the wound.

Highly Absorbent

- Highly absorbent wound pad controls bleeding and absorbs wound exudates

Low adherent

- Special low adherent coating minimizes wound trauma associated with changing dressings

Individually packed and sterile

- Dressings are available individually on demand in an easy tear wrapper and sterile for wound cleanliness.



Ease of application

- The wound pad is attached to a wrapping leader of gauze bandage making the dressing simple to apply, yet comfortable and secure for all general wound dressing.





More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes