



AEROWIPE Povidone Iodine Swabs 60 x 33mm Box/100

\$14.72

In stock
SKU AWP16100

 Add to Wish List

 Add to Compare

Povidone Iodine Swabs (2 ply)



Details

Size: 60 x 33mm **Pricing Unit:** box/100 **Carton QTY:** 100 Our Aerowipe Povidone Iodine swabs are primarily used for large area skin prep prior to injection or surgery.

- 2 Ply Wipe for better skin area coverage
- 10% Povidone Iodine for a gentle antiseptic action



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes