





BLIZZARD Silver Baby Wrap

\$32.08

In stock
SKU BPS-23

 Add to Wish List

 Add to Compare

Reflexcell Baby Wrap - Silver



Details

Size: 65 x 24cm **Pricing Unit:** unit/1 **Carton QTY:** 100 Hypothermia treatment is a critical part of first aid when dealing with trauma induced shock and where the body has insufficient ability to maintain body temperature.

Protect newborns from developing hypothermia immediately after birth with the Blizzard Baby Wrap™. This innovative 3-layer insulating pouch, made of Reflexcell material, delivers Tog 7 level of warmth that is both waterproof and windproof. A newborn may experience difficulty in regulating his/her core body temperature relative to the environmental temperature. The Blizzard Baby Wrap™ Reflexcell technology traps warm air, the metalized surface reflects heat back to the newborn, elasticated triple layer pouch draws material close to the body - reducing cold spaces & convection heat loss, thereby assisting in maintenance of a stable core body temperature. The Blizzard Baby Wrap's resealable double-sided tape closure easily allows paramedics to monitor the newborn while safeguarding the core body temperature. The Blizzard Baby Wrap™ is an important component of an EMS neonatal kit.

FEATURES

- 3-layer pouch
- Reflexcell material
- Double-sided tape closure
- Waterproof
- 7 TOG warmth
- Available color: Orange



- Flat-packed



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes