



CARDIACT Standard AED Cabinet 45 x 35.5 x 14.5cm

\$187.91

In stock
SKU CC-25

 [Add to Wish List](#)

 [Add to Compare](#)


Standard Wall Cabinet for PAD-350P/500P

Details

Size: 47 x 37 x 14.5cm

Pricing Unit: unit/1

Carton QTY: 1

This indoor metal wall cabinet comes with Perspex front and complete with wall brackets. With its ARC approved AED signage, this cabinet ensures that your defibrillator is highly visible, increasing staff awareness and minimising lost time in time of an emergency. This attractive cabinet will also protect your AED from dust and grime. cc25. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes