



CARDIACT Connect Monitored AED Cabinet

\$2,442.82

In stock
SKU CC-300

 [Add to Wish List](#)

 [Add to Compare](#)

Outdoor Wall Cabinet with Monitoring for any AED



Details

Size: -

Pricing Unit: unit/1

Carton QTY: 1



Outdoor Wall Cabinet with Monitoring for any AED



[dfliip id="11144" type="thumb"][/dfliip]



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes