



LAERDAL Adjustable Stifneck Collar - Paediatric (GST FREE)

\$63.04

In stock
SKU CCL-SSC

 Add to Wish List

 Add to Compare

Laerdal Stifneck Collar Select - Adjustable



Details

Size: Paediatric **Pricing Unit:** unit/1 **Carton QTY:** 1 Laerdal Paediatric Stifneck Collar Select - Adjustable Also available in [Adult](#) version.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes