




ChoiceMMed Nylon Carry Case for Pulse Oximeter




\$9.82

In stock
SKU CNC

 Add to Wish List

 Add to Compare

Nylon Carry Case for Pulse Oximeter



Details

Size: Fits C63 and C22 **Pricing Unit:** unit/1 **Carton QTY:** 1 The AERODIAGNOSTIC range offers a selection of commonly used diagnostic supplies in the healthcare industry.

Soft Design

Designed and built with a soft nylon finish.

Durable

Built soft but with a very thick and durable nylon and will keep your oximeter safe

Compatibility

Compatible and designed for most oximeters

Quality

Built with great stitching and a strong nylon material.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes