



## Plastic Screening Kit 36 x 25 x 8.5cm

**\$224.35**

In stock  
**SKU** CPSPK03

 Add to Wish List

 Add to Compare

Plastic Screening Kit 36 x 25 x 8.5cm



### Details

To be used at the screening station prior to staff entry. Contains all the essential products needed for setup and clean-up from conducting the daily scanning procedure. This also contains easy to follow instructions as to how screening is to be conducted.



**KIT CONTENTS:**

<a href="#">ADT10</a>	AERO DIAGNOSTIC Clinical Infrared Forehead Thermometer	1 each
<a href="#">ACG500</a>	AEROCLEANSE Anti-Bacterial Hand Gel 500ml	2 Bottle
<a href="#">AWB7000</a>	AEROWIPE Alcohol Wipes 155 x 140mm	1 tub
<a href="#">AGNPF02</a>	AEROGLOVE Nitrile Glove Powder Free x 2 (Large)	2 5 bag
<a href="#">CM01</a>	AEROMASK Civilian 3Ply Mask W/Ear loops	2 5 each
<a href="#">SG002</a>	AEROSHIELD Safety Spectacles Vortex Clear	2 each
<a href="#">BHB2538</a>	AEROWASTE Clinical Waste Bag 250x300mm - 4L Press 4L Press Seal, 30um	5 bag
<a href="#">CUSTLEAF</a>	Leaflet/Instructions	1 each



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes