



## LIFEPAK CR2 Fully-Automatic Defibrillator with Wi-Fi

**\$3,311.69**

In stock  
**SKU** CR-2-A

 Add to Wish List

 Add to Compare

Fully-Automatic Defibrillator AED WiFi

### Details

**Size:** - **Pricing Unit:** unit/1 **Carton QTY:** 1

- 8-year warranty
- Wireless connectivity
- Fastest time to 1st shock
- CPR Insight
- Durable
- 4 years of standby battery life

Also available as a [Semi-Automatic](#) AED. Sudden cardiac arrest (SCA) can happen to anyone, anywhere. Immediate treatment is vital. A victim's chance of survival dramatically decreases for every minute without treatment. That's why public access defibrillators are so important: they put lifesaving technology where it can do the most good. So when an emergency happens, you should have nothing less than the best. In a cardiac emergency, the LIFEPAK CR2 defibrillator connects everything and everyone involved in the response reducing delays in treatment and improving care across the lifesaving system.


### Who will be saved next?

When you choose the LIFEPAK CR2 defibrillator, you're choosing a connected solution that empowers bystanders to help save a life all while backed by innovative technology that allows continuous CPR and streamlines communication across



the system of care.

With a design influenced by user feedback, it's no surprise the CR2 is rated highest in overall user confidence by AED users. The groundbreaking LIFEPAK CR2 is at the heart of a complete AED response system.

 **Layered design** Layered design with easy-to-follow bold graphics. Both trained and untrained AED users clearly know how to begin. **QUIK-STEP™ electrodes** Peel directly off the base for faster placement. **cprINSIGHT analysis technology** Analyses for shockable rhythm during chest compressions with no need to pause. **Metronome and CPR coaching** Sets an effective pace and audibly guides users. **ClearVoice™ technology** Detects background noise and adjusts volume accordingly for clear instruction. **Fully automatic** Available in fully or semi-automatic models. **Highest available escalating energy** Up to 360J for more effective shocks as needed. **Bilingual** Toggle between two pre-set languages when using the device. **Child mode** Child mode delivers lower energy levels appropriate for young children without having to change electrodes. **LIFEPAK TOUGH™** IP55 rating for challenging environments. **8-year warranty** Backed by an 8-year warranty. **LIFELINKcentral AED Program Manager** Monitor AED programs by tracking AED status, sending patient data to emergency responders and hospitals, detecting AED locations and other tools. [Register your Product Warranty Here.](#)



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes