



## TRAUMASIM Arm Haemorrhage Control Stabilising Bracket

**\$248.05**

In stock  
**SKU** HCS01A

 Add to Wish List

 Add to Compare

TraumaSim Haemorrhage Control Simulator Arm

### Details

This bracket system can be used to secure the TraumaSim Haemorrhage Control Lifecast Arm (HCS01) arm to a bench or table. 2x additional drill holes are available should you wish to permanently secure the bracket and/or arm to a sturdy structure additional bolts & washers are required for this, not supplied.

Includes:

Sturdy & durable plastic bracket

2 x 30mm hex-head bolts

2 x cast iron g-clamps



[Check out the HCS01 TraumaSim Haemorrhage Control Simulator Arm product here](#)



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes