



TRAUMASIM Laceration Haemostatic Clotting Trainer (Square)

\$743.81

In stock
SKU HCT04

 [Add to Wish List](#)

 [Add to Compare](#)

TraumaSim Haemostatic Clotting Trainer Basic Square
Laceration

Details

This medically accurate part-task trainer is made from medical-grade silicone. Improve trainee skills and confidence under pressure using the Haemostatic Clotting Trainer.

The Haemostatic Clotting Trainer range allows users to pack realistic bleeding wounds to practice the vital skill of haemostatic dressing application. Life-cast silicone skin looks and feels real adding realism to improve trainee skill and confidence. These units have the ability to simulate blood flowing from the wound using a pump bottle via tubing. The blood flowing from the wound is stopped when the wound is correctly dressed.

Task trainers allow students to train repeatedly in wound packing and haemostatic dressing techniques via a sturdy, realistic silicone base with tubing for bleeding capability. This is a Laceration Wound with the blood source controlled via a syringe. It is low-cost but critical training aid for any professional delivering life saving first aid training.

Australian designed and manufactured, this simulator is a necessity for any Trainers looking for suitable practical training aids and is guaranteed to add realism to your course delivery!



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes