



BRAYDEN CPR Manikin (Without Lights)

\$516.74

In stock
SKU IM13-S

 Add to Wish List

 Add to Compare

Brayden CPR Manikin (no illuminators)



Details

Size: Adult

Pricing Unit: unit/1

Carton QTY: 1

The Brayden CPR manikin is a realistic training manikin with chest rise and oral & nasal passages, full head tilt and chin lift for opening airway.



Also available in the Brayden [Pulse](#) and [Pro](#).

For our dedicated Brayden Manikin Site, [Click HERE](#).



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes