




Disposable Clear Plastic Kidney Dish 700mL (230mm)

\$1.69

In stock
SKU KD230

 Add to Wish List

 Add to Compare
Plastic Kidney Dish

Details

The AEROSUPPLIES stationary range covers a number of common stationary requirements applicable to the First Aid Industry.

Disposable Clear Kidney Dish

Kidney dishes are a widely used product in surgical procedures for holding small instruments and more importantly for rinsing with sterile water. The dishes are also used to contain samples for laboratory analysis.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes