





HEARTSINE Wire Samaritan AED Wall Bracket 14 x 13 x 13cm



\$172.25

In stock
SKU PAD-CAB-02

 Add to Wish List

 Add to Compare

PAD Wire Wall Bracket



Details

Robust wire wall bracket to suit HeartSine Samaritan AEDs (350P/360P/500P)

Product Features:

- Functional
- Low-cost
- Easy install



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes