



Peak Flow Meter Adult

\$37.04

In stock
SKU PFM01A

 Add to Wish List

 Add to Compare

Peak Flow Meter (Child)



Details

Size: Adult **Pricing Unit:** unit/1 **Carton QTY:** 30 The AERODIAGNOSTIC range offers a selection of commonly used diagnostic supplies in the healthcare industry.

Peak Flow Meters

These devices are required for monitoring patients suffering from respiratory disease or discomfort, to measure oxygen flow rate. **The AeroDiagnostic Peak flow meter is also available in a [pediatric version](#).**



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes