



RAPIDSTOP Small Bleed Control Kit - Tactical 8 x 18 x 6cm

\$138.59

In stock
SKU RSK400


 Add to Wish List

 Add to Compare

RAPIDSTOP Small Bleed Control Kit - Tactical 8 x 18 x 6cm



Details

The RapidStop Tactical Bleeding Control kit is designed for quick and effective control of mass haemorrhaging. Including the popular RapidStop Tourniquet, this haemostatic trauma solution offers military grade ruggedness, all age suitability, superior single-handed application and the fastest average occlusion time. Achieving full occlusion in the shortest time possible is critical. This kit is designed to provide treatment options with efficiency and effectiveness in mind. The Tactical Bleeding Control Kit is compact in size and effective in all environments, whether used in Military, First Response or Civilian instances. 

KIT CONTENTS:

DESCRIPTION	QUANTITY
RAPIDSTOP Tourniquet Orange Strap	1
FIRSTCARE Civilian Trauma & Hemorrhage Control Bandage 10 cm x 17 cm	1
AEROINSTRUMENT Universal Shears 19 cm	1
AEROGLOVE Large Nitrile Powder-Free Glove/ Pair	1
AEROSUPPLIES Marker Pen	1
5.11 TACTICALMED KIT Pouch	1
Dimensions: 20 x 10 x 9cm	



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes