







## AERORESCUE Alloy Scoop Stretcher

\$693.31

In stock  
**SKU** STR-06

 Add to Wish List  
 Add to Compare  
Alloy Scoop Stretcher



### Details

**Size:** 190 x 43 x 6cm **Pricing Unit:** unit/1 **Carton QTY:** 1 AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

### STR-06

**Weight limit:** 160kg Scoop stretcher designed for gently scooping up a patient while maintaining the found position minimising the risk of further injury. Extendable head rest for 170cm-190cm patient size. Folds down to 172cm x 43cm x 6cm



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes