



AEROWASTE Vomit Bag Dispenser

\$66.96

In stock
SKU VBD001

 [Add to Wish List](#)

 [Add to Compare](#)



The AEROWASTE Vomit Bag Dispenser is a wall-mount holder for organised storage of AEROWASTE Bio-Waste Vomit Bags. Designed for clinics, ambulances, schools, aged care facilities and workplace first aid stations - keeps bags accessible, hygienic and visible.



Details

What is the AEROWASTE Vomit Bag Dispenser?

The AEROWASTE Vomit Bag Dispenser is a wall-mount storage holder designed for organised, accessible storage of AEROWASTE Bio-Waste Vomit Bags. Used in GP clinics, urgent care, ambulances, schools, aged care, and workplace first aid stations. Keeps bags visible, hygienic and easy to grab during emergencies.

Where to install

- GP consult rooms and treatment bays.
- Urgent care and triage areas.
- Ambulance patient compartments.
- School sick bays and first aid rooms.
- Aged care nursing stations and residential rooms.
- Workplace first aid stations.
- Public-facing reception desks at high-volume clinics.



Key features

- **Wall-mount design** - hangs from any flat surface with screws or strong adhesive.
- **Holds 25-30 folded AEROWASTE bags** - refill from the 50-pack box.
- **Open-front design** - bags accessible in seconds during an emergency.
- **Durable plastic construction** - cleans easily with disinfectant wipes.
- Compatible with AEROWASTE 1500ml Bio-Waste Vomit Bags.

How to mount

1. Choose a location at adult chest height in clinics, or easy-reach in vehicles.
2. Mark the screw holes on the wall with a pencil.
3. Drill pilot holes and insert wall plugs (supplied).
4. Screw the dispenser firmly to the wall.
5. Load 25-30 folded AEROWASTE bags from a fresh 50-pack box.
6. Pair with a biohazard waste bin within arm's reach.

Related products

- [AEROWASTE 50-Pack refills](#) - bulk box for ongoing supply.
- [AEROWASTE 4-Pack](#) - smaller supply for home/vehicle.
- [First aid kits](#) for complete workplace stocking.

Frequently asked questions

How is the AEROWASTE Vomit Bag Dispenser mounted?

The AEROWASTE Dispenser mounts to any wall, cabinet or vehicle interior with screws (supplied) or strong adhesive strips (sold separately). Position at adult chest height in clinics or in a clear easy-reach location in ambulances and first aid stations. The dispenser stays open for fast access during emergencies.

How many AEROWASTE bags fit in the dispenser?

The dispenser holds approximately 25-30 folded AEROWASTE Bio-Waste Vomit Bags (1500ml). Refill from the box of 50 when stock falls below 10 bags to ensure availability during high-demand periods such as flu season or gastro outbreaks.



Where should an AEROWASTE Dispenser be installed?

Install at every patient-facing point: GP consult rooms, urgent care triage, ambulance patient compartment, school sick bay, aged care nursing stations, and workplace first aid stations. Keep one dispenser per first aid room or vehicle. Always pair with a biohazard waste bin nearby.

Is the dispenser compatible with other vomit bag brands?

The AEROWASTE Dispenser is designed for the AEROWASTE 1500ml bag dimensions. Other brands of similar size may fit but the dispenser is optimised for AEROWASTE products. For best results pair with the AEROWASTE 50-pack refills.

Sources: [Safe Work Australia — First aid in the workplace](#), healthdirect.gov.au.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes