



## QV Gentle Conditioner - 500g

**\$18.41**

In stock  
Only **499** left  
**SKU** 11728

 [Add to Wish List](#)

 [Add to Compare](#)

QV Hair Gentle Conditioner is designed to care for those that experience scalp sensitivity. The unscented conditioner's pH-balanced formulation is dermatologically tested and gently infuses moisture to the hair.

### Details

QV Hair Gentle Conditioner is designed to care for those that experience scalp sensitivity. The unscented conditioner's pH-balanced formulation is dermatologically tested and gently infuses moisture to the hair.

### Directions:

After cleansing with QV Hair Gentle Shampoo, massage QV Hair Gentle Conditioner evenly through the hair, leave for 1 minute, then rinse well. For more intensive conditioning, leave your hair conditioner in for at least 5 minutes before rinsing.

### Benefits:

- Free from fragrance.
- Free from soap and colour.



- pH balanced.
- Ideal for sensitive skin and dry scalps.
- Vegan-friendly - contains no animal-derived ingredients.

## Caution:

Avoid eyes. If product enters eyes flush with water. Store below 30°C.

## Ingredients:

Aqua (Water), Glycerin, Cetearyl Alcohol, Panthenol, Behentrimonium Methosulfate, 1,2-Hexanediol, Caprylyl Glycol, Titanium Dioxide, Citric Acid.