



SES Micro Stretcher with Bag - Orange

\$5,158.27

In stock
Only **100** left
SKU
KITMSORANGE

 [Add to Wish List](#)

 [Add to Compare](#)

In an emergency, fast retrieval of a casualty from a difficult and challenging environment is often life-critical. With the collapsible Micro-Stretcher, it is possible to attend to the injured quickly, safely and effectively in most terrain, be it desert, mountain or battlefield.



Details

In an emergency, fast retrieval of a casualty from a difficult and challenging environment is often life-critical. With the collapsible Micro-Stretcher, it is possible to attend to the injured quickly, safely and effectively in most terrain, be it desert, mountain or battlefield. Assembly from the small carry case is possible in just a few seconds, offering peace of mind for the team and ease of transportation almost anywhere.

Lightweight materials, cutting edge design

The Micro-Stretcher construction comprises of carbon fibre rods with strengthening inserts encased in a hard-wearing netting frame – that contributes to its unique weight, strength and rigidity. The stretcher provides ample support for the casualty with spinal and head support and over-body, chin or head restraining straps. It is tested to withstand up to 311kg* in weight (up to 3.5 times the average UK male). The web mesh material also allows ease of cleaning and is available in a choice of colours to suit its application – military or emergency.

Compact & lightweight

The Micro-Stretcher is a lightweight, compact yet rigid stretcher developed for military and emergency rescue applications over rough terrain or where there may be limited vehicle access. It is ideal for storing in land, marine or aviation vehicles,



and can be utilised as part of a back-pack configuration for ground personnel, allowing them to be hands-free.

Benefits:

- Allows quick & effective assembly in < 20 seconds
- Effective protection & safety - Offering a removable head support with associated strap, 4 colour-coded restraining straps, groin /over body straps & spinal support
- Rapid & easy to assemble - No training necessary
- Adaptable rescue method - Suitable for both horizontal 4-6 person lift & vertical winch by helicopter/ alternative mechanical means
- Ease of securing the casualty - Colour-coded restraining straps and connecting buckles ensures ease of safeguarding in an emergency
- Easy to dismount - Folds into 6 equal sections, allowing ease of storage in its portable case Compact & lightweight portable
- carry case - Protective bag /back-pack or vehicle storage, weighing just 5.1kg in its bag, 4.6kg unpacked
- Strong and durable rigid stretcher - Able to transport a load of 311kg (equivalent to 47 stone)*
- Easy to clean - Mesh fabric & frame can withstand health care cleaning products
- Colour to suit application - standard orange for emergency, military green