



FLUIDSHIELD 3 N95 Particulate Respirator And Surgical Mask - BOX/35

As low as **\$59.87**

In stock

[Only 21 left](#)

Product
availability

Product Name	Qty
FLUIDSHIELD 3 N95 Particulate Respirator And Surgical Mask - BOX/35 - Regular	21.0000
SKU	
FLUIDSHIELD-3- N95-Mask	



Regul
ar

 Add to Wish List

 Add to Compare

Top 5 Reasons to Commit to HALYARD* N95 Respirators

HALYARD* FLUIDSHIELD* N95 Particulate Filter and Surgical Masks

1. 100% of HALYARD* FLUIDSHIELD* N95 Respirators have been cleared by the FDA as a medical device—unlike many N95s or KN95s that have received a temporary Emergency Use Authorization.
2. All HALYARD* N95s are Made in the Americas, unlike suppliers who rely on supply from Asia.
3. All HALYARD* FLUIDSHIELD* Level 3 Respirators deliver Level 3 splash and spray protection — the highest-recognized level of protection.
4. Our duckbill breathing chamber is more than twice as large as the leading competitive surgical N95, and exceeds NIOSH standards for breathability.



5. A balance of proper fit and comfort is critical for health care workers. Only
6. HALYARD* N95s include all of these features:
 1. Two sizes—small and regular
 2. Malleable, adjustable nose wires
 3. Strong elastic straps that are securely bonded—not stapled—to the mask

Details

FLUIDSHIELD* N95 Particulate Filter Respirator And Surgical Mask

- NIOSH-approved
- 5-layer construction
- Soft inner lining
- Protecting healthcare staff since 1997
- 100% natural rubber latex-free
- Available with safety seal
- Designed to provide an effective facial fit
- Easy storage, folds flat without crushing
- 100% of HALYARD FLUIDSHIELD N95 Respirators have been cleared by the FDA as a medical device—unlike many N95s or KN95s that have received a temporary Emergency Use Authorization for the COVID-19 pandemic
- All HALYARD N95's are Made in the Americas, unlike suppliers who rely on supply from Asia
- Our duckbill breathing chamber is more than twice as large as the leading competitive surgical N95, and exceeds NIOSH standards for breathability
- Straps are bonded, not stapled to the mask - Malleable, adjustable nose wire for customized fit

Sold as one Unit BOX/35

SPECIFICATIONS

Brand	FLUIDSHIELD*
UNSPSC Classification	42131606
Sterile	No
Color	Orange
Type	Respirator
Level of Protection	Level 3
Name	FLUIDSHIELD* Ported N95 Particulate Filter Respirator and Surgical Mask
Type	Respirator
FLUIDSHIELD*	Yes
SO SOFT* Lining	Yes



NIOSH	Yes
Fog Resistance	N/A
Fluid Resistant	Yes
ASTM Level of Protection	3
BFE (3)	≥ 99%
STRIKE THROUGH ASTM (80,120,160)	160
PFE (.1)	≥ 99%
Items Per Box	35
Boxes Per Case	6
Items Per Case	210