



Reusable Hot and Cold Pack by Trafalgar

\$15.60

In stock
Only **999** left
SKU 101456

 Add to Wish List

 Add to Compare

Trafalgar's Hot and Cold Packs can be used as both an ice pack and heat pack. The soft pack easily contours to most parts of the body providing a therapeutic relief. Each pack contains 2x reusable hot and cold packs and towelling sleeve.

- Reusable and Microwavable
- Cold Compress & Hot Compress
- Bag Size: 135 x 275mm (WxH)
- Comes with Towelling Sleeve

Hot Compress

You can either heat it up in hot water or the microwave for your convenience. Suitable for **muscular aches, menstrual pain, body stiffness, rheumatism** and **strains**.

Cold Compress

Simply put it in the freezer for at least 1 hour and it's ready to use. Suitable for **headaches, migraines, bruising, swelling, sprains** and **tooth aches**

Details

Reusable Hot and Cold Pack

Reusable and multi-use, great to keep in your home

Trafalgar's Hot and Cold Packs can be used as both an ice pack and heat pack.

The soft pack easily contours to most parts of the body providing a therapeutic relief.



Hot Compress

You can either heat it up in hot water or the microwave for your convenience.

Suitable for **muscular aches, menstrual pain, body stiffness, rheumatism** and **strains**.

Cold Compress

Simply put it in the freezer for at least 1 hour and it's ready to use.

Suitable for **headaches, migraines, bruising, swelling, sprains** and **tooth aches**

Dimensions:

Bag Size: 135 x 275mm (WxH)

Contents:

2x reusable hot and cold packs and towelling sleeve.