



Blue Plastic Adhesive Strips by Trafalgar

\$5.48

In stock
Only **987** left
SKU 101455

 Add to Wish List

 Add to Compare

Ideal for Cuts, Scratches & Abrasions | Suitable for Food Handling

Trafalgar's Blue Plastic Adhesive Strips offers comfort and flexibility for wound care. Sterile, non-stick wound pad and ventilation in the strip to allow the skin to breathe while the wound is healing, suitable for everyday general use.

- Detectable – Visually and electromagnetically
- Water Resistant Adhesive – Seals wound from outside contamination.
- Non-Stick Wound Pad – Protective cushion for the wound.
- Island Dressing – Seals all around, protecting the wound.
- Extra Wide – 25mm wide for extra coverage.
- Sterile and Suitable for Food Handling.

Plastic Strips come in a box of 50

Details

Description:

Trafalgar's Blue Plastic Adhesive Strips offers comfort and flexibility for wound care. Sterile, non-stick wound pad and ventilation which allows the skin to breathe while the wound is healing. The bright blue adhesive makes them easily detectable which is perfect for manufacturing and restaurant industries. Also suitable for everyday general use.

- Ideal for Cuts, Scratches & Abrasions & Suitable for Food Handling
- Detectable – Visually and electromagnetically
- Water Resistant Adhesive – Seals wound from outside contamination.
- Non-Stick Wound Pad – Protective cushion for the wound.
- Island Dressing – Seals all around, protecting the wound.
- Extra Wide – 25mm wide for extra coverage.



- Sterile and Suitable for Food Handling.

Benefits:

- Ventilated – Allows wound to breathe while healing.
- Non-Stick Wound Pad – Protective cushion for the wound.
- Conforming – Flexibility for comfort.
- Island Dressing – Seals all around, protecting the wound.
- Sterile – Individually wrapped for your protection.

Contents:

Plastic Strips come in a box of 50